



eCLIP: Antibody Immunoprecipitation Validation Protocol

Required materials:

1. Cell pellet (non-crosslinked is fine) of 5×10^6 cells
2. Antibody (5 ug per sample plus appropriate amount for western)
3. Magnetic separator (ThermoFisher DynaMag-2 or equivalent)
4. Sonicator (BioRuptor, Qsonica Q800R2 or equivalent)
5. Western blotting equipment
6. Thermomixer (or equivalent; for 65°C sample denaturation)
7. Dark room or ECL Western Imager

Buffer	Provided by	Long-term Storage
eCLIP Lysis Buffer	Eclipse BioInnovations	-20°C
High-Salt Buffer (HSB)	Eclipse BioInnovations	-20°C
25× No-Salt Buffer (NoS Buffer)	Eclipse BioInnovations	-20°C
2× Western Buffer	Eclipse BioInnovations	-80°C

Reagents

Dynabeads M-280 sheep anti-rabbit (for rabbit primary antibodies)*

*for western blot with TruBlot secondary	10 mg/ml	LifeTech	11204D
(or Dynabeads Protein G)	30 mg/ml	LifeTech	37002D

Standard Western blotting reagents (or equivalent):
PVDF membrane

Day 1

Prepare eCLIP Lysis mix

- Pre-chill **eCLIP Lysis Buffer**
- Per sample (5 million cells): add **2.75 µl 200× Protease Inhibitor Cocktail III** to **500 uL eCLIP Lysis Buffer**, mix

Couple antibody to magnetic beads

- **Beads and antibodies:**
 - Use **50 µl beads** per sample
 - rabbit antibodies: use sheep anti-rabbit beads
 - mouse antibodies: use sheep anti-mouse beads
 - others – use Protein G or Protein A beads as recommended
 - Use **5 µg antibody** per sample
- **Prepare beads:**
 - Magnetically separate beads, remove supernatant
 - Wash beads 2x in 500 µl cold **eCLIP Lysis Buffer**
 - Resuspend beads in 100 µl cold **eCLIP Lysis Buffer**



- **Bind antibody:**
 - Add antibody (5 µg) to 50 µl washed beads
 - Rotate, room temp, 45 min

Lyse cells

- Retrieve cell pellets from -80°C freezer
- Immediately add 500µL cold **eCLIP Lysis mix** to each pellet of 5×10^6 cells, pipette to resuspend
- Lyse 5 mins on ice

Sonicate & clarify lysate (while ab+bead binding):

- Sonicate in Bioruptor at 'low' setting, 4°C, 5 min, 30sec on / 30 sec off
- Centrifuge 15,000g, 4°C, 3min
- Transfer supernatant to a new tube

Capture RBP-RNA complexes on beads

- Wash antibody beads 2x in 500 µl cold **eCLIP Lysis Buffer**
- Remove 20 µL (4%) of Sample as backup inputs for western; store at -20°C
 - ✓ This is a backup input sample in case the **input sample** is not properly taken in Day 2
- Add remainder to washed antibody beads
- Rotate 4 degC, 2 h or overnight (in cold room)

Prepare buffers for Day 2

- Dilute 25× **NoS Buffer** according to **Table 1** in a 50mL conical tube. **Mix 1× NoS Buffer and store at 4°C overnight.**

Table 1. 1× NoS Buffer preparation

Component	Volume (mL)
Molecular Biology Grade Water	48
25× NoS Buffer	2
Total volume	50

- Thaw High-Salt Buffer (HSB) at 4°C overnight.



Transfer to membranes

- **Prepare transfer:**

Note: the below is for standard (wet) transfer. We also commonly use rapid transfer methods (such as ThermoFisher iBlot2 system) here as well; if using these systems, follow protocols provided by manufacturer.

- (Have pre-prepared COLD (4 deg) transfer buffer with methanol: 1× NuPAGE transfer buffer, 10% methanol)
- Prepare PVDF membrane(s): pre-flash 10 s in methanol, move to transfer buffer with methanol
- Wet sponges and Whatman papers in transfer buffer with methanol
- Assemble transfer stacks, from bottom to top (black side of stack holder on bottom):
1× sponge – 2× Whatman paper – gel – membrane – 2× Whatman paper – 1× sponge

Transfer:

- overnight 30V (preferred) OR
- 2 hr 200 mA (if doing this, only hook up one transfer box per power supply)

Day 3

Develop COLD membrane

- Block in 5% milk in TBST, room temp, 30 min
- Probe with primary antibody: 0.2 µg/ml (1:5000 for a 1 mg/ml Bethyl antibody stock; can be different per antibody) in 5% milk in TBST, room temp, 1 hr.
- Wash 3x with TBST, 5 min
- Probe with secondary antibody: 1:4000 Rabbit TrueBlot HRP in in 5% milk in TBST, room temp, 1–3 hours
- (Note: if western fails or signal is low, 1:1000 gives higher signal)
- Wash 3× with TBST, 5 min
- Mix equal volumes of ECL Buffer A + Buffer B (or 40:1 of ECL Plus Buffer 1 to Buffer 2), add to membrane and incubate (mix/rotate) for 1-5 min. (1ml final volume per membrane)
- Develop 30 sec & 5 min, then judge signal (15 min maximum; if 15 sec is still too bright, expose two films)