

# 5' End-Seq™

Reveal 5' UTR Landscape and Transcription Start Site Usage Genome-wide

## HIGHLIGHTS

### Define 5' UTR Landscape

Genome-wide detection of known and novel transcription start sites.

### Single Nucleotide End Resolution

Transcription start sites are detected with single nucleotide resolution.

### Identify Alternative Promoters

5' End-Seq detects alternative promoters for all endogenous genes.


## Introduction

The transcription start site (TSS) is the location of transcription initiation at the 5' end of a gene sequence. Knowledge of the exact position of a TSS of an RNA molecule is crucial for the identification of regulatory regions immediately flanking it. 5' End-Seq facilitates active transcript end detection, enrichment and mapping, and identifies known and novel transcription start sites at single nucleotide resolution in only 2 days.

## Define 5' UTR Landscape

End-Seq can facilitate precise UTR identification for RNA therapeutics targeting and can be used as a tool for identifying UTR biomarkers associated with disease. 5' End-Seq can quantify the relative usage of TSS across samples and indicates transcript isoform presence in annotated transcriptomes. Defining the 5' UTR with End-Seq can help predict binding factor motifs more reliably.

## 5' End-Seq Workflow



**Figure 1. 5' End-Seq Workflow.** Total RNA is enriched for polyA-positive mRNA fragments. mRNA is fragmented and uncapped fragments are removed enzymatically. The final library will contain a sense-strand mRNA end fragment, where the 5' end of the read will begin at the 5' end of the mRNA.


## Specifications

<b>Input Sample</b>	Total RNA	>3 ug*
	RNA Concentration	>0.1 ug/ul
	RIN	>7
<b>Sequencing Recommendations</b>	Instrument	Illumina
	Sample Depth	10-20M reads
	Run Parameters	SE100


\*End-Seq service requires submission of 6ug of total RNA

## 5' End Enrichment

5' End-Seq enriches the sample for reads in the 5' UTR 4-fold higher than the coding region. Greater depth of coverage around TSS allows for more precise detection of known and novel TSS.



**Figure 2.** Comparison of read coverage across all hg38 gencode v35 transcripts. RNA-Seq sample shows even coverage across transcripts, while 5' End-Seq library shows enrichment at the 5' end.




**Figure 3.** Distribution of reads in RNA-Seq and 5' End-Seq. Majority of reads in RNA-Seq samples are found in the coding region (CDS) of a transcript, while less than 10% of reads are found in the 5' UTR. Reads in the 5' UTR constitute almost 20% of total 5' End-Seq reads.

## Ordering information

More information about 5' End-Seq services online at [www.eclipsebio.com](http://www.eclipsebio.com) or contact us at [info@eclipsebio.com](mailto:info@eclipsebio.com).

## Genome-Wide End Calling


5' End-Seq sequencing data can be used to call TSS genome-wide at single nucleotide resolution. RNA-Seq samples do not allow for calling TSS, as the read coverage does not change dramatically around each TSS. Easily calling TSS at single nucleotide resolution genome-wide may facilitate biomarker and therapeutic target identification.



**Figure 4.** Example of transcription start site calling in RNA-Seq sample and 5' End-Seq sample for the gene TRABD. Transcription start sites are called confidently in 5' End-Seq sample, and not defined in the RNA-Seq sample.

## Precise TSS Discovery

5' End-Seq is not only able to identify transcription start sites (TSS) with single nucleotide resolution, but the TSS end calls from 5' End-Seq data are near annotated TSS. RNA-Seq provides data spread across the entire gene but unable to call the 5' ends with the same precision as 5' End-Seq.



**Figure 5.** Distance from 5' End-Seq end call to the closest annotated TSS. End calls are found within bases or at annotated TSS. Applying an end call approach to RNA-Seq data identifies less precise 5' ends in fewer genes.