

ECLIP: PREPARATION AND UV CROSSLINKING OF ADHERENT CELLS

Required materials:

1. UV crosslinker with 254-nm wavelength UV bulbs (UVP CL-1000 Ultraviolet Crosslinker or equivalent)
2. Liquid nitrogen or dry ice (sufficient to submerge tubes, in appropriate container)
3. 1x DPBS (*Corning cat# 21-031-CV or equivalent*)
4. Trypan blue stain (*ThermoFisher Scientific, cat# 15250-061* or other equivalent live cell counting assay)
5. Standard cell counting system (hemocytometer or automated cell counter)

Cell viability validation (prior to crosslinking):

- a. Use Trypan blue stain (*ThermoFisher Scientific, cat# 15250-061* or other equivalent live cell counting assay to assay cell viability)
- b. Cell viability should be > 95% to ensure intact RNA

Wash cells:

- a. Aspirate spent media
- b. Wash the plate gently with chilled 1x DPBS (15mL for a 15cm plate)
- c. Aspirate media
- d. Add enough chilled 1x DPBS to just cover the plate (5mL for a 15cm plate)
*** Note:** If all plates are at equal cell density, one plate can be sacrificed for counting – this plate would be dissociated (with trypsin, accutase, or equivalent) and cell number (per plate) counted at this stage. **This is recommended for cell types that require chemical dissociation enzymes** to dissociate and be properly counted, as this is not recommended post-crosslinking

UV crosslinking:

- e. Place the tissue culture plate on leveled ice or a cooling block pre-chilled to 4°C
- f. Place the above (plate plus ice or cooling block) into the UV cross-linker
 - Notes: **Ensure the plate is leveled**
 - **Remove** tissue culture plate lid **before** cross-linking
- g. Cross-link at 254-nm UV with an energy setting of 400 mJoules/cm²

- **Note:** this is a setting of **4000** on many cross-linkers which display values in 0.1 mJoules/cm²
- h. While keeping the cells on ice, use a cell scraper (*Corning cat# CLS3010-10EA*) to scrape the plate
 - i. Transfer the cells to a 15mL conical tube
 - j. Wash plate once with 5mL of chilled 1x DPBS and add to the same 15mL tube
 - k. Centrifuge the 15mL conical tube at 250 x g for 5 minutes at 4 °C
 - l. Aspirate and discard supernatant
 - m. Resuspend in the desired amount of chilled 1x DPBS for flash freezing
 - Typically, 10×10⁶ cells per mL
 - n. Transfer desired amount into 1.5mL Eppendorf Safe-Lock Tubes (or equivalent)
 - Typically, 1mL of 10×10⁶ cells per mL
 - o. Spin down at 250 x g for 5 minutes at 4 °C
 - p. Aspirate the supernatant and freeze by submerging the lobind-tubes completely in liquid nitrogen or on dry ice
 - q. After frozen (at least 30 seconds), remove from the liquid nitrogen (or dry ice) and store at -80 °C