

PREPARATION AND UV CROSSLINKING OF ADHERENT CELLS

Required materials:

1. UV crosslinker with 254-nm wavelength UV bulbs (UVP CL-1000 Ultraviolet Crosslinker or equivalent)
2. Liquid nitrogen or dry ice (sufficient to submerge tubes, in appropriate container)
3. 1x DPBS (*Corning cat# 21-031-CV or equivalent*)
4. Trypan blue stain (*Thermo Fisher Scientific, cat# 15250-061* or other equivalent live cell counting assay)
5. Standard cell counting system (hemocytometer or automated cell counter)

Cell viability validation (prior to crosslinking):

- a) Use Trypan blue stain (*Thermo Fisher Scientific, cat# 15250-061* or other equivalent live cell counting assay to assay cell viability)
- b) Cell viability should be > 95% to ensure intact RNA

Wash cells:

- a) Aspirate spent media
- b) Wash the plate gently with chilled 1x DPBS (15mL for a 15cm plate)
- c) Aspirate media
- d) Add enough chilled 1x DPBS to just cover the plate (5mL for a 15cm plate)

***Note:** If all plates are at equal cell density, one plate can be sacrificed for counting – this plate would be dissociated (with trypsin, Accutase, or equivalent) and cell number (per plate) counted at this stage. **This is recommended for cell types that require chemical dissociation enzymes** to dissociate and be properly counted, as this is not recommended post-crosslinking

UV crosslinking:

- a) Place the tissue culture plate on leveled ice or a cooling block pre-chilled to 4 °C
- b) Place the above (plate plus ice or cooling block) into the UV cross-linker

*** Note:** Ensure the plate is level. **Remove** tissue culture plate lid before cross-linking.

- c) Cross-link at 254-nm UV with an energy setting of 400 mJoules/cm²
- *Note:** this is a setting of 4000 on many crosslinkers which display values in 0.1 mJoules/cm²
- d) While keeping the cells on ice, use a cell scraper (Corning cat# CLS3010-10EA) to scrape the plate
- e) Transfer the cells to a 15mL conical tube
- f) Wash plate once with 5mL of chilled 1x DPBS and add to the same 15mL tube
- g) Centrifuge the 15mL conical tube at 250 x g for 5 minutes at 4°C
- h) Aspirate and discard supernatant
- i) Resuspend in the desired amount of chilled 1x DPBS for flash freezing
- Typically 20×10⁶ cells per mL
- j) Transfer desired amount into 1.5mL Eppendorf Safe-Lock Tubes (or equivalent)
- Typically 1mL of 20×10⁶ cells per mL
- k) Spin down at 250 x g for 5 minutes at 4°C
- l) Aspirate the supernatant and freeze by submerging the LoBind tubes completely in liquid nitrogen or on dry ice
- m) After frozen (at least 30 seconds), remove from the liquid nitrogen (or dry ice) and store at -80°C