

PREPARATION AND UV CROSSLINKING OF SUSPENSION CELLS

Required materials:

1. UV crosslinker with 254-nm wavelength UV bulbs (UVP CL-1000 Ultraviolet Crosslinker or equivalent)
2. Liquid nitrogen or dry ice (sufficient to submerge tubes, in appropriate container)
3. 1x DPBS (*Corning cat# 21-031-CV or equivalent*)
4. Trypan blue stain (*Thermo Fisher Scientific, cat# 15250-061 or other equivalent live cell counting assay*)
5. Standard cell counting system (hemocytometer or automated cell counter)

Cell viability validation (prior to crosslinking):

- a) Use Trypan blue stain (*Thermo Fisher Scientific, cat# 15250-061 or other equivalent live cell counting assay*) to assay cell viability)
- b) Cell viability should be > 95% to ensure intact RNA

Preparation of suspension cells:

- a) Pool all cells per biosample (if multiple plates)
- b) Transfer cells with media to 50mL conical tube(s)
- c) Centrifuge at 250 x g for 5 minutes at 4 °C

Wash cells:

- a) Aspirate spent media
- b) Resuspend the pellet(s) in 25mL of chilled 1x DPBS
- c) Count cell concentration (either with automated cell counter or hemocytometer)
- d) Spin down remaining sample in 50mL conical tube(s) at 250 x g for 5 minutes at 4 °C
- e) Aspirate supernatant
- f) Resuspend cells to no more than 12×10^6 cells per mL

UV crosslinking:

- a) Aliquot a minimum of 10×10^6 and a maximum of 60×10^6 cells (re-suspend in 1x DPBS) in at least 5mL total volume to a standard 15cm tissue culture grade plate

* **Note:** Ensure the cells are evenly dispersed and the plate is fully covered (5mL should be a sufficient volume)

- b) Place the tissue culture plate on leveled ice or a cooling block pre-chilled to 4 °C
- c) Place the above (plate plus ice or cooling block) into the UV crosslinker

* **Note:** Ensure the plate is level. **Remove** tissue culture plate lid before cross-linking.

- d) Cross-link at 254-nm UV with an energy setting of 400 mJoules/cm²

***Note:** this is a setting of 4000 on many crosslinkers which display values in 0.1 mJoules/cm²

- e) After crosslinking is completed, transfer the cells to a 15mL conical tube
- f) Wash plate once with 5mL of chilled 1x DPBS and add to the same 15mL tube
- g) Centrifuge the 15mL conical tube at 250 x g for 5 minutes at 4 °C
- h) Aspirate and discard supernatant
- i) Resuspend in the desired amount of chilled 1x DPBS for flash freezing
 - Typically 20×10⁶ cells per mL
- j) Transfer desired amount into 1.5mL Eppendorf Safe-Lock Tubes (or equivalent)
 - Typically 1mL of 20×10⁶ cells per mL
- k) Spin down at 250 x g for 5 minutes at 4 °C
- l) Aspirate the supernatant and freeze by submerging the LoBind tubes completely in liquid nitrogen or on dry ice
- m) After frozen (at least 30 seconds), remove from the liquid nitrogen (or dry ice) and store at -80 °C